

Features

- Multi-modality: photon-counting, 2D imaging, 3D time-of-flight ranging, TCSPC (time-correlated single-photon counting)
- Image dimension: 32x32 (1024) pixels
- In-pixel counter: 6 bit (photon-counting)
- In-pixel TDC: 10 bit (photon-timing)
- Max frame rate: 100,000 fps (burst) and 10,000 fps (continuous)
- Timing resolution: 312 ps – 0.9 ns
- Full scale range: 320 ns – 0.92 μ s
- Hardware interface: USB 2.0
- Software interface: Matlab

Applications

- Single-photon counting, timing and imaging
- Medium- and Long-distance LIDAR
- Medical imaging, FLIM, FRET
- Time-resolved Spectroscopy
- Diffuse Optical Tomography
- Looking around the corner
- Advanced 3D Gesture user interface
- Safety and Security



Fig. 1: SPAD camera for 2D imaging, 3D ranging and TCSPC photon-counting.




Fig. 2: Array chip containing the 32x32 pixels with 1024 SPAD detectors and 1024 TDCs.

General Description

The TCSPC (Time-Correlated Single-Photon Counting) camera (Fig. 1) is based on 32x32 smart-pixels, each containing a SPAD (Single-Photon Avalanche Diode) detector and a Time-to-Digital Converter (TDC). The chip (Fig. 2) is fabricated in a CMOS technology [1].

Each one of the 1024 pixels provides single-photon detection sensitivity (Fig. 3) and directly measures the time delay between a synch pulse (e.g. triggering a laser shining the sample under test) and the detection of a return single photon. The time delay can be either the photon time-of-flight in 3D-ranging/LIDAR applications, or the fluorescence photon emission in FLIM/FRET imaging, or the time-resolved waveform of very fast optical signals. The 1024 fully independent pixels operate in parallel, with no multiplexing during either detection or TDC conversion, in a global shutter mode. The pixel noise is the best-in-class (Fig. 4) among CMOS SPAD detectors.

In 2D/intensity mode, at every frame the SPAD array provides 1024 words of 6 bits each (i.e. single-frame photon-counting dynamics up to 64 photons). In 3D depth-resolved/timing mode, at every frame the SPAD array provides 1024 words of 10 bits each (i.e. single-frame time-bin down to 312 ps and 320 ns full-scale range).

The frame rate tops at 100 kfps (i.e. a minimum integration frame-time of 10 μ s) in burst mode, and 10 kfps in continuous mode.

For more information please contact:

Prof. Franco Zappa

email: franco.zappa@polimi.it

Phone: +39-02-2399-6149

Dr. Alberto Tosi

email: alberto.tosi@polimi.it

Phone: +39-02-2399-6174

www.everyphtoncounts.com

Electrical Characteristics with $V_{EX} = 5 V$

Parameter	Notes	Min	Typ	Max	Units
Resolution	number of pixels	32 x 32			pixels
SPAD diameter	diameter of the SPAD into each pixel	30			μm
Pixel pitch	size of each pixel	150			μm
Pixel fill-factor	with no microlens	3.14			%
Photon Detection Efficiency (PDE) of SPAD detectors	@ 420 nm	50			%
	@ 300 nm and @ 650 nm	20			%
	@ 800 nm	5			%
Dark count rate (DCR)	with $V_{EX} = 5 V$	100			cps
Afterpulsing	with 20 ns hold-off time	2.6			%
Optical cross-talk	between nearest pixels	0.01			%
Photon-counting dynamics	number of bits of the in-pixel counter	6			bit
Photon-timing dynamics	number of bits of the in-pixel TDC	10			bit
On-board clock (f_{clock})	adjustable through the GUI interface	70	200		MHz
TDC least significant bit (LSB) ^[1]	with $f_{\text{clock}} = 200 \text{ MHz}$		312		ps
	with $f_{\text{clock}} = 160 \text{ MHz}$		390		ps
	with $f_{\text{clock}} = 70 \text{ MHz}$		900		ps
TDC full-scale range (FSR) ^[1]	with $f_{\text{clock}} = 200 \text{ MHz}$		320		ns
	with $f_{\text{clock}} = 160 \text{ MHz}$		400		ns
	with $f_{\text{clock}} = 70 \text{ MHz}$		920		ns
Full-scale depth range	max scene depth, around the target distance (even kilometers from the camera, depending on the optical setup)	48	60	138	m
TDC single-shot precision (FWHM)	with $f_{\text{clock}} = 160 \text{ MHz}$ (all 32x32 TDC working in parallel)		660		ps
Frame-rate	Burst mode (limited by FPGA read-out)	no limit		100,000	fps
	Continuous mode (limited by USB 2.0 link)	no limit		10,000	fps
Camera power supply			+5		V
Lens connector		12 mm – F/1.4 C-Mount			

^[1] LSB and FSR are adjustable through the GUI by changing the on-board clock frequency (f_{clock}).
 $1 \text{ LSB} = 1/(f_{\text{clock}} \cdot 16)$ $n = 10 \text{ bit}$ $\text{FSR} = 2^n \cdot \text{LSB} = 1024 \cdot \text{LSB}$




Fig. 3: Photon Detection Efficiency at $V_{EX} = 5 V$ of the CMOS SPADs.




Fig. 4: Dark Count Rate distribution at $V_{EX} = 5 V$ and at room temperature of the CMOS SPADs.

[1] F. Villa, D. Bronzi, Y. Zou, C. Scarcella, G. Boso, S. Tisa, A. Tosi, F. Zappa, D. Durini, S. Weyers, U. Paschen, and W. Brockherde, "CMOS SPADs with up to 500 μm diameter and 55% detection efficiency at 420 nm," *J. Mod. Opt.*, pp. 1–14, Jan. 2014.




Fig. 5: Pixel block diagram with 30 μm SPAD and 10 bit TDC.




Fig. 6: Single-shot precision with 312 ps bin (left), and 390 ps bin (right). The FWHM precision is 600 ps and 660 ps, respectively.




Fig. 7: Centroid (left) and FWHM (right) of the arrival time histogram (about 40,000 conversions) in all the pixels of the array (bin = 312 ps).





Fig. 8: Centroid (left) and FWHM (right) of arrival time histograms of each array pixel with low (red) and high (blue) optical power illuminating the 32x32 SPAD + TDC array.




Fig. 9: SPAD + TDC camera Graphical User Interface for 2D imaging and 3D ranging.




Fig. 10: 3D photon-timing ranging image acquired by the 32x32 SPAD camera at $\lambda = 750 \text{ nm}$ (90 mW uniformly distributed over the acquired scene) in 1 ms (dwell time). Left image at 6 m, right image at 8 m.




Fig. 11: Target at 8 m acquired with a standard camera (left), with the SPAD camera in 2D photon-counting mode (center) and in 3D photon-timing mode, with the reconstruction merging 2D and 3D data (right). Scene was scanned with 25 steps (i.e. reconstructed image has 160 x 160 pixels) for a total of 250 ms dwell time.




Fig. 12: Target at 8 m acquired with a standard camera (left) and reconstruction by merging 2D and 3D data (right) acquired by the SPAD camera. Scene was scanned with 50 steps (i.e. reconstructed image has 160 x 320 pixels) for a total of 250 ms dwell time.




Fig. 13: Frames from a 2D movie of a circular neon lamp, acquired at 50,000 fps, showing the fast gas discharge propagation through the tube.